
Beipackzettel (BPZ) für Hüftprothesen (DE) Instruction Leaflet Hip Total Prosthesis (EN) Foglietto Illustrativo per Protesi Totale di Anca (IT)

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EN

1. General Security Information

Products from Atesos medical SA (Atesos) should exclusively be used by surgeons, who got an equivalent instruction and have experi­ences in implanting Total Hip Prostheses. Atesos recommends attending product specific training and / or demonstrating surgeries per­formed in reference hospitals. Before using Atesos products the operating surgeon should study this package insert, the product labelling and the information provided in the surgical technique. Current documentation about these products is available from the sales force or online at www.atesos.ch.. The manufacturer is not liable for complica­tions which arise from incorrect indication or in­correct surgical technique, incorrect selection of implant components, incorrect use of the im­plant-related instruments or non compliance to aseptic rules.

Contraindications

Before the operation the surgeon is bound to inform the patient about the course of the sur­gical procedure, about the involved risks, con­traindications and side effects as well as the necessary measures of precaution (see also Chapter 3. product specific information). The patient should be advised, how his own behav­iour influences the security and longevity of the joint replacement. The surgeon should docu­ment all information provided to the patient.

Implant Passport

The surgeon has to provide the patient with an implant passport. The patient should preserve and carry this document. It serves to provide information for other attending physicians and as a document to pass safety locks using metal detectors. In case of loss the patient should ask his physician for replacement of the docu­ment.

Hypersensitivity of the patient towards im­plant materials

Hypersensitive reactions or allergic reactions caused by implant materials, in particular metal ions have to be considered. The labelling on each packaging specifies for each implant the kind of material used. The attending physician should consider in his preoperative planning the potential risks arising from hypersensitive reactions caused by implant materials.

Reutilisation or Modifications

Implants are intended for single use only and should only be used in original condition with intact packaging. No modifications are al­lowed; exempt these are specified in the surgi­cal technique.

Cleaning and sterilisation, Resterilisation

Atesos implants are sterilised by Gamma or X-ray radiation (min 25 kGy) or by Ethylene- Oxide. Implants which are delivered as sterile must not be resterilised by the customer. Im­plants and Instruments, which are delivered as non-sterile have to be cleaned and sterilised and can be resterilised if needed.

Implants and instruments designated for rest­erilisation, should not get into contact with ag­gressive or deleterious substances. For Steri­lisation instruments and implants to be sterilised should be placed in a suitable containment. In­struments must not be sterilised in the bags they are packed for delivery. The refurbish­ment and steam sterilisation has to be per­formed according to the valid standards (ISO 17664). Further information for disinfection, cleaning and sterilisation of instruments see brochure "Refurbishment of reusable Instru­ments" which is provided by Atesos as part of the product documentation. In case of sus­pected TSE contamination the residence time should be 18 minutes (134°C).

Storage & Transport

The implants should be stored and transported dry at room temperature or at situational tem­perature in the original package non opened and protected against direct sunlight. Beyond that there are no further particular require­ments for transport, storage or handling. After the expiry- date indicated on the label they must not be used.

2. Product specific information

The implant components are intended for a Partial or Total Hip Joint replacement. They are intended for primary or revision interventions according to surgical technique.

Indication:

- Primary and secondary arthritic degenera­tion of the hip- joint,
- non-vascular necrosis of the femoral head
- rheumatoid arthritis provided sufficient bone quality
- dysplastic cox arthritis up to 145° CCD an­gle
- posttraumatic alterations provided sufficient bone quality
- distinct Bone diseases, such as Achondro­plasia
- For revision products: revision interventions (for more detailed information see corre­sponding surgical technique).
- promising joint-preservative reconstruction

Contraindications

- Larger deformations and defects of the fe­mur or the Acetabulum.
- Circulatory disorders in the neck of the fe­mur in avascular neck roses (MRI scan)
- Radiological manifest osteoporosis or Oste­omalacy, locally or general.
- Metabolic disorders that can affect bone for­mation
- progressive local or general tumor develop­ment,
- x-ray aggravated bone substance, (with the exception of perioperative radiation for ossi­fication prophylaxis)
- acute infection of the joint or the environ­ment, passed or still eminent infection dis­eases,
- Severe muscle, nerve or vascular diseases that endanger the extremity
- Material hyper-sensitivity
- Allergy to metals (titanium)
- Modified medullary canal geometry is un­suitable for the stem geometry
- Infected skin lesions on the operating leg

Relative Contraindications

- Infected skin lesions on the contralateral leg
- Sustained and still threatening infectious disease with possible joint manifestation
- pregnancy
- Diabetes mellitus with an HbA1c value of> 7mmol / l

Risk factors for the success of the hip joint replacement

- systemic diseases and metabolic insuffi­ciencies
- Infections in Anamnesis
- downfalls in Anamnesis
- drug, nicotine-, alcohol-, pharmaceutical-misuse
- obesity, (Body Mass Index >30)
- epilepsy
- serious sport, strong physical strain
- Thrombosis or pulmonary embolism during surgery
- postoperative morphological alterations of the load bearing structures or of the involved tissue may lead to implant failure caused by: loosening, deformation, crack formation or fracture of the component (or cement) or bone early wear of the implants

Potential adverse reactions, potential complications

- dislocation, subluxation, restricted range of motion, undesirable shortening or lengthen­ing of the leg
- infection
- pain
- vein- thrombosis and pulmonary embolism
- functional disorder of the heart-lung system, of the lungs or nerve system (e.g. due to fat embolism)
- hematoma, delayed wound healing
- bone- fractures due to high or unequal bone loading or weakening of the bone substance
- loosening or migration of the implant
- wear of the articulating surfaces may induce osteolysis
- fracture of the implant, bone or cement.

These risk factors and complications may re­quire the Reoperation of the artificial joint.

Information about use and combination of Atesos products

The physician should verify that the selected components are compatible to each other and he has to consider the limitations in combina­tion of materials and devices. All stems and cup products of the Atesos range can be combined with products officially dis­tributed by Atesos. Combinations with other third-party products are at the user's risk, must be checked and, if necessary, authorized For Ceramic-Ceramic articulation only inlays and ball-heads of the same manufacturer or ap­proved combinations can be combined. Ball-heads made of stainless steel (ISO 5832-1) are not approved for use with stems made of unalloyed Titanium (ISO 5832-2) or stems made of any CoCr alloy. Femoral heads may only be combined with femoral stems from Atesos, if the combination is approved by Atesos. The cone connection stem – ball head and inlay - shell may only be mounted if the cone surfaces are undamaged and clean. Standard ceramic heads are only approved if the stem- cone is new (unused). If in case of a revision the stem is not going to be replaced, only ceramic ball-heads with titanium

sleeves or metal ball-heads are allowed as a replacement. In case of a broken ceramic ball-head before, solely revision ceramic ball heads with titanium sleeves are accepted to be used. In case of revisions damaged ball heads or inlays have to be exchanged. In the case of re­vision surgery, care must be taken to ensure that the implant surfaces of the implants remaining in the patient are not damaged. Damage in­creases the risk of breakage.

The stems have been tested according to ISO 7206-4 and -6 in combination with ball heads up to the lengthiest neck length XXL It is in the responsibility of the surgeon, if a larger neck length is being used.

The surgeon must be familiar with the fixation concept of the implants, the basics for plan­ning and the surgical technique. The indica­tion and selection of the implant sizes and models and the implant positioning influence early and late complications (early implant mi­gration, dislocations, loosening and increased wear).

Combination of Instruments

For the implantation the implant related surgi­cal instruments of Atesos have to be used. Ex­empted are instruments that are generally used for joint surgery or that are described in the surgical technique. In case of unauthorised combination of products of other manufac­turers Atesos denies all liability

3. Surgical Procedure

The implantation has to be performed respect­ing the attached surgical technique. For the success of the surgical procedure it is crucial that the surgeon is familiar with the recom­mended surgical technique.

Preoperative Planning

The preoperative planning enables the sur­geon to select the appropriate implant com­ponents, and to decide upon the combination within the range. The procedure should be planned on the base of medical findings includ­ing the imaging based diagnostics. In case that another but the planned size will be indicated or the planned implant cannot be used, other implant sizes or models have to be available. Failures in preoperative planning and insuffi­ciencies in imaging can cause improper im­plant type or implant size selection or cause improper implant positioning. For further infor­mation concerning preoperative planning and the use of digital x-ray templates consult the surgical technique or contact Atesos or the re­sponsible sales force.

Preparation of the surgical procedure

Instruments:

The instruments are delivered non sterile and have to be cleaned and sterilised before use. Before use the proper function of the instru­ments should be verified. The use of damaged instrumnts may cause complications or early failure of the joint replacement.

Implants:

When unpacking the implants the information on the label and on the implant should be checked (Ref, Model, Size and Sterility Date). During unpacking the principles of asepsis have to be respected. For tracking purposes the LOT- no of the used implant should be docu­mented in the patient's record. For that pur­pose additional stickers are included in the packaging.

Any contact with other items that may damage or contaminate the surface has to be avoided. Before implantation the components have to be checked for integrity. Protection covers should be removed just before use.

Implant components must not be im­planted if

- damaged or scratched
 - inappropriately treated or if used not ac­cording to the surgical technique
 - implanted before
 - packaging or packaging seal is damaged
- Please send damaged implants to the Atesos sales force or to the manufacturer.

Intraoperative use

For implants which are cemented the operating physician has to respect the instructions of the manufacturer of the cement.

Beipackzettel (BPZ) für Hüftprothesen (DE)
Instruction Leaflet Hip Total Prosthesis (EN)
Foglietto Illustrativo per Protesi Totale di Anca (IT)

Summary of Safety and Clinical Performance available on <https://ec.europa.eu/tools/eudamed>
502.C001-D7 Ed. 11/22 / CE1250

The initial intraoperative fixation decides on the primary stability of the implant and is an essential factor for the long term success of the joint replacement.

The following factors may initiate or encourage implant early loosening or other complications:

- improper preparation of the bone seating
- selection of inappropriate implant size
- incorrect positioning of the implants
- improper primary stability
- fractures or other bone damages

Before closure the operative field should be cleaned of cement, bone and other particles.





Postoperative Treatment

For the postoperative treatment, accepted procedures have to be applied and the requirements formulated in the surgical technique have to be respected. The patient has to be informed about the limits during rehabilitation and later on of the joint replacement in terms of loading, range of motion and sports activities. Periodic follow-up examinations of the postoperative course are to be carried out to suit the individual and general risk situation. In case of appropriate evidence, examinations at shorter intervals may be necessary.

Interaction of implants with X-rays and magnetic fields


Patients, who have a hip joint implanted from Atesos, can be diagnostically treated without negative consequences, either with radiation based (conventional X-rays /CT) or magnetic based (MRI) techniques. Depending on the type of equipment, the quality of the imaging will be differently influenced.

4. Symbol description

	Manufacturer
	European Authorized Representative
	Respect instruction manual
	Single use only
	Do not re-sterilise
	Do not use product after expiry date (Year – Month – Day)
	Do not use if package or seal is broken
	Sterilised in the final packaging through radiation
	Sterilised in the final packaging through Ethylene- Oxide
	Single sterile barrier system with outer protective packaging
	Double sterile barrier system
	Non-cemented use
	Reference number / Article number
	Batch number of the manufacturer
	Medical device
	Unique Device Identifier

References: If you need additional information about the use of our product read the surgical technique of Atesos and / or refer to your sales force. Disposal information: Before disposing of instruments and implant components that have been in contact with patients, disinfect and sterilize them (see also ISO 12891 Surgical Implants - Removal and Analysis of Surgical Implants). The disposal takes place according to material-specific aspects, according to the specifications of the institution.

Obligation for incident report: All users are obliged to report serious incidents in accordance with the product to the responsible authorities.

	
Atesos medical AG Schachentallee 29 CH-5000 Aarau Tel: +41 (0)62 823 15 15 Fax: +41 (0)62 823 26 94 https://atesos.ch	hestomed GmbH Am Eichenwald 13 D-09356 St. Egidien Tel: +49 37204 638-0 https://www.hestomed.de
	

	
Atesos medical AG Schachentallee 29 CH-5000 Aarau Tel: +41 (0)62 823 15 15 Fax: +41 (0)62 823 26 94 https://atesos.ch	hestomed GmbH Am Eichenwald 13 D-09356 St. Egidien Tel: +49 37204 638-0 https://www.hestomed.de
	

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